

Lightmouse H®

USE INSTRUCTIONS
Lightmouse Home



What is *Lightmouse Home*?

Lightmouse Home emits an infrared athermic light and it is a device dedicated to the self aesthetic treatment for both women and men. It is apt to treat skin blemishes of the face, neck, decollete', breasts and it is also valid for the man in case of problems of scarce or partial stiffening of the penis.

Light, source of life

Light is a source of life: both the plant and animal kingdoms get strength and wellness from light. *Lightmouse Home* is based on this simple concept: use light to bio-stimulate the vital functions while respecting the integrity and health of the body. The light emitted by *Lightmouse Home* is infrared at 905 nm (thus not visible by the human eye) and athermic (= cold).



Applications

important:

Read carefully the "GENERAL REMARKS, INSTALLATION AND GUARANTEE" sheet enclosed in the packaging.

In no case may the *Lightmouse Home* device and its applications be intended as a cure or diagnosis practices. This device is apt only to improve the aesthetics of the part treated. It is not a medical device.

note:

- Use the disposable bag for each application with *Lightmouse Home*. Once finished, you can buy them on or else you can use a plastic film for food preservation.
- Wash thoroughly the area before treating it with *Lightmouse Home* as any substance (lotions, oils, etc) becomes a barrier for the penetration of the light.
- In order to appreciate the results, treat the area daily for at least 10-15 days. As for the duration of each application, follow the instructions of this manual.

Note about the "cold" infrared light

Although there is an energy supply, Infrared at 905-30 nm can be considered cold as it does not increase the temperature of biological matter. This feature guarantees that the device is safe and there are no side effects, for ex. burns. The feeling of warmth that may be felt at the end of a treatment is due to the improvement of the blood circulation and lymphatic circulation in the tissues.

References

In the scientific world, Low Level Laser (or LED) Therapy, or the so called LLLT, is a method well known to treat and solve several problems. Light at 905 nm combined with the proper device is considered and used as anti-inflammatory therapy by traditional medicine (physiotherapy).

Scientific studies and bibliographies report the real benefits that such methods produces in the body.

The technology

The emitters *Lightmouse Home* are equipped with 16 LEDs (Light Emitted Diode). This technology ensures long durability and constant emission and the end user can be sure they will get the maximum benefit with no risks at all. Besides, the disposable bag is not only a guarantee of hygiene but also a safety guarantee as it functions as insulation between the emitter and the user.